

The Agilent 2100 Bioanalyzer

RNA Integrity Number (RIN)

A standardized approach for RNA integrity assessment

qPCR 2005
Freising
Marc Valer



Acknowledgment

Software

- Andreas Schroeder
- Michael Leiber
- Susanne Stocker
- Thomas Ragg



Application

- Samar Lightfoot
- Ruediger Salowsky
- Marcus Gassmann
- Christine Miller*
- Rainer Wittig⁺

Array Correlation work

- Kohei Ando


We would like extend special thanks to our collaboration partners Ambion Inc. and the German cancer research center as well as to Quantiom Bioinformatics. We would also like to thank all researchers who beta-tested the RIN and provided valuable inputs.

* Johns Hopkins Medical Institute

⁺ German Cancer Research Center (DKFZ)


RNA LabChip kits

Analysis of Total RNA Integrity



Typical first QC step after RNA sample prep prior to microarrays or real-time PCR

High quality total RNA



Partially degraded total RNA

← 2100 bioanalyzer: single lane gel-like image

Current total RNA integrity measures

UV offers no information on size distribution, and does “see” interferences of DNA, solvents and other contaminants. Sophisticated de-convolution is possible to individually quantify the components.

Gels, have a limited linear range and ribosomal band intensity ratio had been used to differentiate good to bad samples. This measure is very subjective, dependent on the

RNA degradation is a gradual process.

Results have to be interpreted by visual inspection.

Overlay of electropherograms only works well for samples with the same concentration.

Instrument dependency in signal height

mRNA QC is addressed with smear analysis tools

2100 expert - C:\...iles\Agilent\2100 bioanalyzer\2100 expert\data\samples\demo\electrophoresis\Demo mRNA Pico.xad (Read-Only)

File Context View Electropherogram Windows Help

Data Overlaid Samples

Electropherogram - Sample 7

Assay Properties Chip Summary Gel Electropherogram Result Flagging Log Book

All Files

- Demo mRNA Pico.xad
 - All Samples
 - Sample 1
 - Sample 2
 - Sample 3
 - Sample 4
 - Sample 5
 - Sample 6
 - Sample 7
 - Sample 8
 - Sample 9
 - Sample 10
 - Sample 11
 - Ladder

Sample 7

	From [nt]	To [nt]	Corr. Area	% of Total	Color
1	200	2,966	193.7	67	Red
2	3,077	12,942	81.1	28	Blue

Local Global

Advanced Expand


- Alignment
- Quantitation
- Sizing
- RNA Fragment
- Smear Analysis
 - Perform Smear Analysis
 - Regions T.
- Baseline calculation
- Filter Settings
- Baseline correction
- Integrator

Default

Results Peak Table Fragment Table Region Table Legend Errors

Auto Export Auto Print Auto Run

Problem Description



The ratio of ribosomal bands is not sufficient to describe RNA integrity!

RNA degradation is a gradual process.

Results have to be interpreted by visual inspection.

Overlay of electropherograms only works well for samples with the same concentration.

Instrument dependency in signal height

The RNA Integrity Number (RIN)

The RIN is a software algorithm that allows classification of RNA integrity


It extracts a number of characteristic features from the bioanalyzer electropherogram

An adaptive learning process is used to “teach” the algorithm about the relative importance of the extracted features (based on many example electropherograms).

The RIN returns a value that is characteristic for the integrity of a specific sample.

Feature Generation


Approach: Training of regression models for peak positions and Intensities based on peak labels



Features
Peak height/position
Areas/Area ratios
S/N ratio
Max, Min values
Waviness of the curve
...


Selection of the best feature combination

RIN Visualization



Anomaly threshold settings

: RNA Integrity Number	
Pre Region Anomaly Thres...	0.5
5S Region Anomaly Thresh...	0.3
Fast Region Anomaly Thre...	0.77
Inter Region Anomaly Thre...	0.21
Precursor Region Anomaly ...	0.3
Post Region Anomaly Thre...	0.29
Baseline Anomaly Threshold	0.56
Ribosomal Ratio Anomaly T...	0.89
Unknown Sample Type Thr...	0.83
Marker Anomaly Threshold	0.88



RIN - Limitations

What the RIN can do:

- **Obtain an assessment of RNA integrity.**
- **Directly compare RNA samples**
- **Ensure repeatability of experiments**


What it CANNOT do:

- **Predict the outcome of an experiment if no prior validation was done**

Results



RIN Application – Assessment of RNA Integrity




Intact RNA: RIN 10

Partially degraded RNA:
RIN 5

Strongly Degraded RNA:
RIN 3

RIN Application – Directly Compare Samples

same sample on different instruments




When testing an identical RNA sample on various instruments, identical RINs are obtained – within narrow limits

**36 samples 3 instruments
CV RIN: 3 %
CV ribosomal ratio: 12 %**



RIN Application – Directly Compare Samples

same sample in different dilutions



When testing an identical RNA sample in various dilutions, identical RINs are obtained – within narrow limits

**108 samples 3 dilutions
CV RIN: 3 %
CV ribosomal ratio: 17 %**

RIN Application – Comparison to alternative integrity measures

When applying the degradometer tool to a random data set, a significant number of samples are “unclassified” (BLACK CODING)


integrity class	% labeled BLACK
10	21
9	30
8	27
7	35
6	20
5	20
4	42
3	44
2	100
1	100

Unclassified samples

The degradometer is based on the ribosomal ratio and the area corresponding to our “fast region”.

The degradometer classifies samples in 4 categories, when RIN data is clustered in blocks the data correlates.

RNA QC in Routine Gene Expression Workflow



Experimental: Correlation of RIN to Microarray

Materials and Methods






- Preparation of Total RNA Degradation
- Start Total RNA Samples
- Design of Experiments

Results

- Labeled cRNA Quality Check by BioAnalyzer
- Scatter Plots : self to self
- Scatter Plots : Intact vs. Degradation
- Reappearance : Compare Plots


Conclusions

Material and Method

Steps	Material & Reagent	Conditions
<u>Preparation of Total RNA</u>	Biochain Human Brain Total RNA	Human Brain Intact Total RNA Brain Partial Degradation Brain Complete Degradation
<u>Labeling & cRNA check</u>	Agilent Low RNA Input Linear & Amplification Kit (5184-3523) 	Total RNA 500 ng start Agilent Standard Protocol
<u>Microarray</u>	Agilent Human 1A (ver. 2) (G4110B) 	Agilent 2100 BioAnalyzer (RNA 6000 Nano LabChip kit) NanoDrop 
<u>Hybridization & Scanning</u>	In Situ Hybridization Kit Pls (3184-3568) SureHyb Chamber (G2534A)  Agilent Scanner (G2565BA)	SSC Wash Ozone Free 

Preparation of Total RNA Degradation

Products from BIOCHAIN Co.
Human Brain Total RNA
Intact




60 °C
8 hours




60 °C
2 hours


Partial Degradation



Complete Degradation



Start Total RNA Samples : Intact, Partial Degradation and Complete Degradation



- Intact
- Partial Degradation (60C 2 hours)
- Complete Degradation (60C 8 hours)

	260/230	260/230	RIN	rRNA Ratio
Intact	2.14	1.94	7.3	1.3
Partial Degradation	2.11	1.85	5.2	0.5
Complete Degradation	2.13	1.78	2.4	0.0

Design of Experiments



Samples : Human Brain Total RNA (from BIOCHAIN Co.)

- Intact
- Partial Degradation
- Complete Degradation


No.	Barcode	Cyanine 3	Cyanine 5	Comment
1	251209734174	Intact	Intact	Self
2	251209735151	Partial Degradation	Partial Degradation	Self
3	251209735054	Complete Degradation	Complete Degradation	Self
4	251209735149	Intact	Partial Degradation	Differential
5	251209735150	Partial Degradation	Intact	Differential
6	251209735053	Intact	Complete Degradation	Differential
7	251209735055	Complete Degradation	Intact	Differential

} Dye Swap


} Dye Swap

Labeled cRNA Quality Check by BioAnalyzer

Cy3 labeled cRNA



Cy5 labeled cRNA




Peak Size : 0.1K (Complete Deg.)


- Ladder
- Intact (RIN = 7.3)
- Partial Degradation (RIN = 5.2)
- Complete Degradation (RIN = 2.4)




Scatter Plots : Self to Self




Intact




Partial Degradation




Complete Degradation




Scatter Plots : Intact vs. Degradation




Intact (self to self)



Intact vs. Partial Degradation



Intact vs. Complete Degradation





Other developments

The principles applied for the RNA integrity number are valid for any kind of curve where differences, even if slight, are highly reproducible. The algorithm is ideal for sample classification where categories can be predefined and are pure.

A new algorithm which is able to automatically select the most informative features has been in the meantime developed. It proved already a good correlation with the expert selected RIN features. It will now be applied for other classification tasks

RIN to DIN? aCGH sample QC amplification

Genomic DNA amplification with phi29 titration series (50ng, 0.05ng, 0.005ng & no input) analyzed with the RNA 6000 Nano kit. The increasing non-specific peaks can be seen in the 30-40 second range and the decreasing specific peaks can be seen migrating in the 50-60 second range.



50ng genomic DNA

0.05ng

No Input

The same samples were analyzed on a 0.8% agarose gel and visualized with EtBr. 600ng of each sample was loaded
Lanes: 1-1kB marker, 2-50ng, 3-5ng, 4-0.5ng, 5-0.05ng, 6-0.005ng, 7-0.0005ng, 8- no input
Can not distinguish between the specific and non-specific amplification peaks.

Summary

- **A software algorithm was presented that allows standardization of RNA integrity assessment**
- **Ribosomal ratio is important parameter for RNA integrity but many other factors have to be taken into account**
- **Results indicate that the software measurements are independent from the instrument variability and can ensure repeatability of results (e.g. RT-PCR or microarray)**
- **Degradation level of starting RNA material seriously influences microarray results.**
- **Labeled cRNA from different levels of RNA degradation results in low Cy Dye incorporation and low yield of cRNA**

Thank you for you attention!!!

Q&A...all along

Live Demo:

- Kit preparation
- Instrument setup
- Sample: RNA Ladder, Ambion Human Liver totalRNA
- Start Run
- See results